


PATENTED TECHNOLOGY

SaeboMAS Mini

PRODUCT MANUAL



Maximizing function while improving proximal strength.



www.Saebo.com

Symbols Used

Symbols used in this Product Manual

SYMBOL	MEANING
	Attention! Read these parts extra carefully for they contain important information regarding safety

Symbols used on the SaeboMAS

SYMBOL	MEANING
	SaeboMAS Mini is a Medical Device
	Declaration of Conformity with EU MDR
	Manufacturer's name and address
	Production date; Country of Manufacture
	Serial Number
	Unique Device Identifier
	Catalogue Number
	Warning; Crushing of hands. Closing mechanical parts can injure hands
	Warning; Moving parts can crush
	Refer to instruction manual/booklet
 WARNING To avoid risk of injury reduce tension to "1" prior to securing or unsecuring this bar to the frame	Warning! Reduce tension to the lowest level prior to securing or unsecuring this bar to the frame, in order to avoid risk of injury
 WARNING Secure only to stable surfaces Avoid plastic or wheeld tables	Warning! Secure the table clamp only to stable surfaces. Avoid plastic tables or tables with wheels

Contents

Symbols Used 2

Introduction 4

 Possible Benefits 4

 Contraindications..... 4

 Precautions 4

Parts List 5

Warnings 5

Accessories List..... 6

 Pronation/Supination Support 6

 Purpose Wrist Support..... 6

 Purpose Elbow Support..... 6

Installation 7

 A. Attach the Table Mount 7

 B. Attach SaeboMAS Mini to Table Mount 8

 C. Height Adjustment 8

 D. Apply Forearm Support..... 9

 E. Apply Elbow Support..... 9

 F. Apply Wrist Support 10

 G. Forearm Support Lock..... 10

 H. Tension Adjustments..... 11

 I. Correct Position 12

 J. Swivel Stop..... 12

Start Using SaeboMAS Mini..... 13

 Troubleshooting Tips 13

Cleaning 14

Maintenance..... 14

Warranty 15

Introduction

The **SaeboMAS Mini** is a lightweight, **zero gravity, dynamic mobile arm support** that challenges and assists the weakened shoulder and elbow during functional tasks. This compact device was specifically designed for clients suffering from neurological and orthopedic conditions resulting in shoulder and elbow weakness. It is the ideal home program solution to improve strength, motor recovery, and independence. Clients will be able to perform functional exercises with greater ease and minimal compensation. Additionally, this personal device allows for enhanced independence for self-care, leisure or occupational tasks like using a computer, eating, drinking, or grooming.

The main feature of the SaeboMAS Mini is the tension system. Based on the functional goals, the adjustable spring based parallelogram provides precise customization for assistance or resistance. Whether the client requires extra assistance for functional tasks or less assistance to further challenge his or her strength during exercise, the proprietary tension system is fully adjustable.

POSSIBLE BENEFITS	CONTRAINDICATIONS
• Increase motor control and strength • Increase range of motion • Improve self-care performance • Minimize over-use injuries and unwanted movement • Provide an opportunity for the patient to perform highly repetitive tasks	The SaeboMAS Mini should NOT be used with individuals who exhibit severe shoulder and elbow pain

PRECAUTIONS

To decrease risk of bodily injury, please be sure to lower the Tension Indication Knob completely until you feel a stop (knob parallel with bottom arrow) prior to opening or closing the SaeboMAS Mini (clipping or unclipping).

The SaeboMAS Mini must be mounted to a sturdy table. The device should not be mounted to unstable surfaces or tables on wheels. Do not mount to plastic tables.

Do not perform hand-to-mouth activities without the Elbow Support. The Elbow Support must be attached. Without the Elbow Support attached, the Forearm Support may migrate distal (i.e., slip downwards) towards the wrist. In some cases, the Forearm Support may slip completely off the arm.

For some individuals, the Forearm Support may migrate distally (i.e., slip downwards) toward the wrist during use. Add non-slip material, such as waffled shelf liner to minimize migration.

Parts List





Figure 1

WARNINGS

SUBJECT	WARNING
Amount of support	 Too much support can result in shoulder pain
Swivel stop	 Be sure to loosen the Swivel Stop Knob prior to use to avoid damage


Accessories List



PRONATION/SUPINATION SUPPORT

SaeboMAS Mini Pronation/Supination Support is for SaeboMAS and SaeboMAS Mini users who want to have more forearm movement. The lightweight design allows for easy pronation and supination movements of the forearm.

Important: SaeboMAS Mini Pronation/Supination Support is only meant to be used with the SaeboMAS and SaeboMAS Mini. Read the SaeboMAS or SaeboMAS Mini Product Manual before using and assembling the SaeboMAS Pronation/Supination Support.



PURPOSE WRIST SUPPORT

The Wrist Support is ideal for patients that exhibit “wrist drop”.

PURPOSE ELBOW SUPPORT

It is highly recommended that the elbow support be used for any activities involving elbow flexion. For example, all activities that incorporate hand to mouth movements (i.e., feeding, grooming) will need the elbow support. This will prevent the Forearm Support from migrating distally (slip downwards) towards the wrist.

Installation

During the installation, please keep the **SaeboMAS Mini** in a closed position (i.e., horizontal bar secured/clipped to the frame) until the setup is complete (see **Figure 3**). Once you are finished installing the device, unclip the horizontal bar from the frame of the SaeboMAS Mini to begin treatment. When you are not using the device, lower the Tension Indication Knob completely until you feel a stop (knob parallel with bottom arrow) (see **Figure 4**).

Important: To decrease the risk of bodily injury, make sure the tension is lowered completely before opening/closing the SaeboMAS Mini (see **Figure 4**).



Figure 3




Figure 4

A. ATTACH THE TABLE MOUNT

1. The Table Mount should be positioned on the involved side approximately 2 feet (60 cm) away from the user (See **Figure 5**).
2. The Table Mount should be secured to a sturdy heavy surface. Ideally, the mounting surface should be 1-3 inches (2.5 - 7.5cm) thick. If you are applying the Table Mount to a counter top, be sure to have a sufficient surface area to properly secure the mount to the counter.
3. Place the Table Mount in the correct position (see **Figure 6**).
4. Tightly rotate the Table Mount Handle until the mount is secured (see **Figure 7**).




Figure 5



Figure 6



Figure 7

Installation (Cont'd.)

B. ATTACH SAEBOMAS MINI TO TABLE MOUNT

1. Now that the Table Mount is secured to the surface, insert the SaeboMAS Mini vertical pole into the opening of the Mount (see Figure 8).
2. Once the appropriate height is determined, secure the vertical pole by rotating the Height Adjustment Knob on the Table Mount in a clockwise direction (see Figure 9).
3. Be sure to keep the device in a closed position (see Figure 3 on page 6).



Figure 8



Figure 9

C. HEIGHT ADJUSTMENT

To increase or decrease the height of the device, rotate the Height Adjustable Knob in a counter-clockwise direction (to loosen) and gently pull up or push down to the appropriate height. Once the appropriate height is determined, secure the vertical pole by rotating the Height Adjustable Knob in a clockwise direction (see Figure 9).

Important: What height should the SaeboMAS Mini be positioned? The Tension Adjustment Handle of the device should be positioned at shoulder height of the user (see Figure 10).




Figure 10

Installation (Cont'd.)

D. APPLY FOREARM SUPPORT

The SaeboMAS Mini includes 2 Forearm Supports (1 Small, 1 Medium/Large)

1. To remove the Forearm Support, loosen the Knobs of the Supports Clamp and remove the Forearm Support.
2. To attach a Forearm Support, place the Forearm Support in the slot of the Supports Clamp (see Figure 11), and slide it to the other end of the Supports Clamp. Re-secure the knobs of the Supports Clamp.
3. Forearm Support position should be approximately 1/3 of the forearm-length to the elbow. (See Figure 12)



Figure 11



Figure 12

Reminder: The Forearm Supports are malleable and can be shaped for an intimate fit.

E. APPLY ELBOW SUPPORT

1. To attach the Elbow Support, loosen the Supports Clamp Knobs. Slide the Elbow Support into the slot on the side as indicated in Figure 13. Re-tighten the Knobs to secure the Elbow Support to the Supports Clamp.
2. Be aware that elbow extension will be limited while using the Elbow Support.

Important: **X** Do not use the elbow support for tasks that require elbow extension.
✓ Make sure that the Forearm Base is still fully through its slot.



Figure 13

Installation (Cont'd.)

F. APPLY WRIST SUPPORT

To attach the Wrist Support, loosen the Supports Clamp Knobs and slide the Wrist Support into the slot on the side as indicated in **Figure 14**. Re-tighten the Knobs to secure the Wrist Support to the Forearm Support.

Important: Make sure that the Forearm Base remains fully through its slot




Figure 14

G. FOREARM SUPPORT LOCK

One of the unique features of the SaeboMAS Mini is the ability to lock and unlock the forearm support as needed. For example, patients that exhibit poor strength at the elbow and forearm may require the forearm support to be in a locked positioned for increase stability and control (i.e., Forearm Support unable to swivel vertically/up and down). Conversely, patients that have fair or good strength at the elbow and forearm may be able to tolerate the forearm support in an unlocked position (i.e., Forearm Support able to swivel vertically/up and down).

1. To lock the Forearm Support, rotate the Knob until it latches in the slots. Rotate the Forearm Base a bit back and forth till it locks in place. (See **Figure 15**).
2. To unlock the Forearm Support, pull the black knob and rotate the knob 90 degrees.




Figure 15

Installation (Cont'd.)

H. TENSION ADJUSTMENTS

The SaeboMAS Mini includes a Tension Scale to identify and track the amount of support needed for the affected arm (see **Figure 16**). The scale ranges from min to max (min=least amount of tension/support; max=greatest amount of tension/support).




Figure 16

1. Tension Adjustments
 - a. Increase Tension/Support: Rotate the Tension Adjustment Handle in a clockwise direction.
 - b. Decrease Tension/Support: Rotate the Tension Adjustment Handle in a counter-clockwise direction.
 - c. To avoid damage, once the desired tension is reached, be sure to move the Tension Adjustment Handle away from the frame.
2. Progressing the Program
 - a. When appropriate, consider decreasing the amount of tension provided to further challenge the affected arm.

Important: How much support should be provided? Each individual will require a customized program based on his or her motor impairments. There should be enough support provided to successfully complete the functional task. At the same time, the support should be limited to challenge the user. Finding a balance between too much support and too little will be important.

Installation (Cont'd.)

I. CORRECT POSITION

The SaeboMAS Mini should be used in a seated position. The user should be positioned side-by-side (approximately 2 feet or 60cm away) from the SaeboMAS Mini. (See Figure 17)




Figure 17

J. SWIVEL STOP

Apply the Swivel Stop to eliminate rotation or turning of the frame when the SaeboMAS Mini is not in use. Be sure to unlock prior to use.

1. **Unlock the Swivel Stop:**
 - a. Pull the black knob and rotate the knob 90 degrees (See Figure 18)
2. **Lock the Swivel Stop:**
 - a. Rotate the black knob till it latches back into closed position (Figure 19)
 - b. Rotate the SaeboMAS Mini till the Swivel Stop is aligned with the yellow hole (Figure 20)
 - c. At this location, the Swivel Stop will automatically lock in place



Figure 18



Figure 19



Figure 20

 **Important:** Be sure to unlock the Swivel Stop Knob prior to use to avoid damage.

Start Using SaeboMAS Mini

To start working with the SaeboMAS Mini, follow the steps below:

1. Determine the right height and position of the SaeboMAS Mini relative to the user
2. Decide which size Forearm Support and which accessories are needed
3. If applicable, mount the accessories correctly
4. Check if the accessories are securely mounted
5. Check if the Tension Indication Knob is at the lowest position
6. Loosen the Swivel Stop
7. Unclip the Horizontal Bar
8. Put the user's forearm in the Forearm Support and apply the Forearm strap
9. Apply the right compensation level

To stop working with the SaeboMAS Mini follow the steps below:

1. Turn the compensation level to minimal
2. Unstrap the arm
3. Clip in the Horizontal bar
4. Tighten the Swivel Stop

TROUBLESHOOTING TIPS

1. If the user exhibits decrease wrist extension strength (i.e., wrist drop) consider applying the Wrist Support for added support.
2. During elbow flexion activities, the elbow may shift/migrate while using the Elbow Stop. For these instances, consider securing the elbow to the Elbow Support by using a Velcro strap.

Cleaning

When the SaeboMAS Mini is successively used by more than one user, parts that have been in contact with the patient shall be disinfected before the successive user uses the SaeboMAS Mini.

When cleaning the outer part of the SaeboMAS Mini, make sure to use a soft, clean cloth dampened with lukewarm water and a non-abrasive liquid household cleaner.

 **Important:** Avoid cleaning the cable and the bearings.

Maintenance

Forearm liners are removable and can be replaced by ordering directly from Saebo. To increase the longevity of the liners, consider applying a stockinette over the client’s arm prior to using the device.

Check the cable integrity regularly. Decreasing cable integrity can be recognized by cable fraying. In case cable fraying is observed, contact your distributor.

 **Important:** The SaeboMAS Mini must, under no circumstances, be opened or dismantled in order to clean or check inner parts. If opening is required, contact your distributor or the manufacturer.

Warranty

Limited Warranty: The **SaeboMAS Mini** was manufactured to meet Saebo’s superior standards. However, regardless of how carefully you use, or how well you care for your mobile arm support, it will eventually begin to show age and wear. This warranty does not cover damage caused by accident, improper care, negligence, normal wear and tear, or natural breakdown of colors and materials over extended time and use. Saebo, Inc. warrants that this product will be free from physical defects for a period of one (1) year from the date of purchase. If your product does not work properly because of a defect in materials or workmanship, Saebo, Inc. (referred to as “the warrantor”) will, for the length of the period, which starts with the date of original purchase (“warranty period”), at its option either (a) repair your product with new parts, or (b) replace it with a new product. The decision to repair or replace will be made by the warrantor. This warranty is given in lieu of all other warranties, written or oral, however expressed. Any and all warranties, other than this one, whether expressed or implied including implied warranties or merchantability and fitness of purpose, are hereby disclaimed.

No Liability for Damages: In no event shall Saebo, Inc., or its suppliers, be liable for any monetary damages whatsoever, including any consequential, incidental or indirect damages. This warranty does not apply if the equipment has been subject to misuse, neglect, accidents, improper installation, or put to use in violation of instructions furnished by Saebo; has been repaired or altered outside a repair facility factory sanctioned by Saebo. In the event this warranty is deemed inapplicable by Saebo, a repair estimate including parts, labor and shipping costs shall be provided to the customer, and shall be paid the customer prior to commencement of repair work.

 **Important:** Any severe personal injury in relation to your SaeboMAS Mini should be mentioned to the manufacturer.



www.Saebo.com



Manufactured and Distributed by:
Saebo, Inc.

2459 Wilkinson Blvd, Ste. 120-B
Charlotte, NC 280208, USA

888-284-5433
sales@saebo.com
www.saebo.com

Rev. 3 | 2026



WARNING

This product can expose you to Diisononyl Phthalate (DINP), which is known to the State of California to cause cancer. For more information go to WWW.P65Warnings.ca.gov

EU REP

European Healthcare & Device Solutions (Ireland) Ltd
Stratton House, Bishopstown Road Cork, T12Y9TC, Ireland
info@europeanddevicesolutions.eu
www.europeanddevicesolutions.co.uk

UK REP

European Device Solutions Ltd
5 Nuthatch Close, Wideopen,
Newcastle Upon Tyne, England, NE13 6NT

Australian Sponsor:
Asumedic Australia Pty Ltd
Unit 3, 3 Basalt Road Pemulwuy NSW
2145 1300 473 4222
performhealth.com.au

GMDN: 48071
BUDIDI/GMN: 081001835SaeboMASminiB7

Attn: Intended user group, both laypersons and healthcare professionals.