


SaeboStretch® | Measurement Form

Dynamic Hand Splint



Indications

- Has mild to moderate tone or soft tissue shortening.
- Can passively achieve neutral wrist extension with composite finger extension with minimal effort.



With fingers together, measure straight across the digits (proximal phalanges) as shown in the photo. Avoid wrapping the measuring tape.

XXS: 45 – 51 mm

XS: 51 – 57 mm

SM: 57 – 70 mm

MD: 70 – 83 mm

LG: 83 – 102 mm

Hand: Left: Right:

Size: XX Small: X Small: Small: Medium: Large:

If undecided about the size, assess finger length.
Longer fingers? Consider a larger size.

Total Length of Splint	
XXS	23cm
XS	25.5cm
SM	29cm
MD	33cm
LG	36cm