

Saebo Hand Treatment Kit



Intensify Your Treatment with Evidence-Based Solutions

The Saebo Hand Treatment Kit is a comprehensive upper extremity program designed for hospitals and clinics that treat neurological clients with limited hand function. The all-in-one kit allows clinicians to intensify their treatment by addressing hand function, sensorimotor recovery, spasticity and soft-tissue management.



SaeboGlove

A low-profile device helping clients with neurological and orthopedic injuries incorporate their hand functionally in therapy and at home.



SaeboGlide

Ideal for exercising the effected upper limb. Handstrap allows for easy use with limited hand activity.



SaeboFlex

A functional, dynamic orthoses designed to allow clients with upper extremity spasticity to grasp and release in therapy and at home.



SaeboStretch

A comfortable, dynamic resting hand splint designed to prevent joint damage, maintain bone alignment, and improve range of motion.



SaeboStim Pro

A portable, two-channel electrical stimulation unit providing four types of stimulation: NMES for strengthening & edema, TENS for pain treatment, SES for sensory recovery and cortical priming, and FES for functional training using the built-in trigger button.



SaeboStim One

A simple wireless electrical stimulation device for shoulder subluxation as well as other parts of the body.



Saebo Mirror Box

A scientifically proven treatment tool used to speed up and improve motor function and sensation and decrease the severity of stroke and other neurological disorders.